

A vertical blue bar on the left side of the page, featuring a sunburst or lens flare effect that radiates from the top left corner towards the center.

K-Lift

Age Management System

K E R S T I N  
F L O R I A N

# Energize. Lift. Transform.

## K-Lift Takes Skin Rejuvenation to a New Level

This advanced system integrates the strength of three modalities of Kinetic energy – energy in motion – which work in harmony with the body’s natural energy to target the underlying causes and factors associated with aging.

Visible Light Energy (Photobiostimulation)

Electron Energy (Impulse MicroCurrent)

Magnetic Energy (Transdermal Delivery)

## Targets all layers of the skin for visible results...

Lifts and Tones Facial Muscles

Increases Collagen Production

Reduces Lines and Smooths Wrinkles

Tightens Facial Contours



## Multi-Phase Rejuvenation

Integrating the energy technologies achieves multilayer skin rejuvenation throughout the three treatment phases:

### Phase I/K-Light

#### Stimulates Energy Points and Warms the Skin

Utilizes a 640 nm Red Spectrum Light that energizes the mitochondria, the power center of the cell to restore cell function to younger levels, accelerates the cell’s repair mechanism and improves the cell detoxification process.

### Phase II/K-Firm

#### Tones and Lifts the Muscles, Stimulates Collagen Production

Impulsed MicroCurrent technology triggers DNA regeneration and cell repair, re-educates and rejuvenates muscle tissue and redefines facial contours. Automatically adjusts to skin resistance.

### Phase III/K-Fusion

#### Transdermal Delivery Penetration of Key Nutrients into the Skin

A cutting-edge technology that infuses essential nutrients into the skin to promote collagen regeneration, skin oxygenation and lymphatic circulation.

## K-Lift Synergistic Technologies

K-Light	K-Firm	K-Fusion
Photobiostimulation + Impulse MicroCurrent	Impulse MicroCurrent	Transdermal Delivery
640 nm Red Spectrum Light	Triggers DNA regeneration and cell repair	Utilizes a Bi-polar balanced wave current to help penetrate topical ingredients into the skin
Increases collagen production up to 5x (NASA)	Rejuvenates muscle tissue	
Restores cell function to younger levels	Tones and lifts muscles	Applies cutting-edge technology to oxygenate the skin by increasing blood flow and improving lymphatic circulation
Accelerates cell repair	Redefines facial contour by stimulating muscle receptors	
Improves cell membrane permeability		
Increases detoxification		

# Treatments

## **Youth Accelerator Facial, 60 minutes**

Transform your skin with this modern age-management treatment. Includes a deep cleansing and exfoliating multi-acid peel for the face and neck. Receive 40 minutes of LED Red Light and Impulse MicroCurrent to help awaken cell function to younger levels, re-educate muscle tone and trigger DNA regeneration and cell repair. Treatment concludes with the infusion of essential nutrients into the deep layers of the skin using cutting-edge Transdermal Delivery technology.

## **Advanced Lifting Facial, 90 minutes**

Take your skin to the next level of transformation with this modern age-management treatment. Includes a deep cleansing and exfoliating multi-acid peel for the face and neck. Receive 60 minutes of LED Red Light and Impulse MicroCurrent to help awaken cell functions to younger levels, accelerate repair mechanism, re-educate muscle tone and trigger DNA regeneration and cell repair. Treatment concludes with the infusion of essential nutrients into the deep layers of the skin using cutting-edge Transdermal Delivery technology.

## **Express Treatments, 20 minutes each**

### Express Lift

This preventive treatment is designed to maintain youthful cell energy, strengthen muscle tone and help delay the aging process. Recommended after a series of 6 Youth Accelerator or Advanced Lifting Facials to maintain the results.

### Express Eyes

This targeted treatment provides immediate lifting and brightening around the eye area as it helps to diminish the look of fine lines and wrinkles.

### Express Lips

Lip lines, lower cheeks and jaw area are toned and smoothed with this specialty treatment designed for the lip area. Fine, feathery lines are smoothed and jaw muscles are lifted to help achieve a more youthful-looking appearance.

### Express Congestion

This treatment is designed to help draw out excess sebum from oily skin while increasing lymphatic drainage to reduce puffiness and enhance circulation.